CHROMATOGRAPHIC CONDITIONS:

		<u>ASSAY</u>				
Column No:	169	Type of Column:	250*4.6.Agilent waters			
Column Temp (°C):	25					
Detection λ (nm):	230	Injection Vol (μL):	10			
Mobile Phase: Compos	ition (% v/v	v) & Ratios				
ACN:BUFFER(1:24)			Flow Rate (mL/min):	1.5		
			Pump Pressure (bars):	200-250bars		
DISSOLUTION						
Column No:	0	Type of Column:	UV			
Column Temp (°C):	0					
Detection λ (nm):	uv	Injection Vol (μL):	UV			
Mobile Phase: Compos	ition (% v/v	z) & Ratios				
•	111011 (70 17) & ratios	Flow Rate (mL/min):	UV		
UV			Pump Pressure (bars):	UV		
			rump rressure (burs).			
						

REFERENCE SUBSTANCES:

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.	AMOXICILLIN TRIHYDRATE	A52-1	99.4
2.			
3.			
4.			
5.			