## **CHROMATOGRAPHIC CONDITIONS:**

		<u>ASSAY</u>		
Column No:	100	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	30	<u></u>		
Detection $\lambda$ (nm):	254	Injection Vol (μL):	20	
Mobile Phase: Compos	ition (% v/v	e) & Ratios	Flow Rate (mL/min):	1.5
1598:400:2 PH 5.9			Pump Pressure (bars):	122-125
			<u> </u>	
		DISSOLUTIO	<u>ON</u>	
Column No:	100	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	30	<u></u>		
Detection $\lambda$ (nm):	254	Injection Vol (μL):	20	
Mobile Phase: Compos <u>H2O:ACN:TFA</u> 1598:200:2 PH 5.9	ition (% v/v	) & Ratios	Flow Rate (mL/min): Pump Pressure (bars):	1.5 122-125

## **REFERENCE SUBSTANCES:**

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			