

CHROMATOGRAPHIC CONDITIONS:

ASSAY

Column No: 176 Type of Column: Hyperclone BDS C18, 15cm
Column Temp (°C): 30
Detection λ (nm): 275 Injection Vol (μ L): 20

Mobile Phase: Composition (% v/v) & Ratios

Buffer pH 6.0: ACN(70:30) Flow Rate (mL/min): 1.0
 Pump Pressure (bars): 150-153

DISSOLUTION

Column No: 176 Type of Column: Hyperclone BDS C18, 15cm
Column Temp (°C): 1
Detection λ (nm): 1 Injection Vol (μ L): 10

Mobile Phase: Composition (% v/v) & Ratios

Buffer pH 6.0: ACN(70:30) Flow Rate (mL/min): 1.5
 Pump Pressure (bars): 120-130

REFERENCE SUBSTANCES:

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			