

CHROMATOGRAPHIC CONDITIONS:**ASSAY**

| | | | |
|-------------------|-----|---------------------|--------------------|
| Column No: | 196 | Type of Column: | Xterra RP 18, 25cm |
| Column Temp (°C): | 30 | | |
| Detection λ (nm): | 272 | Injection Vol (μL): | 20 |

Mobile Phase: Composition (% v/v) & Ratios

| | | |
|---------------------------------------|-----------------------|-----|
| Water:MeoH: ACN: GAA(380:820:100:10) | Flow Rate (mL/min): | 1.0 |
| | Pump Pressure (bars): | 180 |
| | | |

DISSOLUTION

| | | | |
|-------------------|-----|---------------------|--------------------|
| Column No: | 196 | Type of Column: | Xterra RP 18, 25cm |
| Column Temp (°C): | 30 | | |
| Detection λ (nm): | 272 | Injection Vol (μL): | 20 |

Mobile Phase: Composition (% v/v) & Ratios

| | | |
|---------------------------------------|-----------------------|-----|
| Water:MeoH: ACN: GAA(380:820:100:10) | Flow Rate (mL/min): | 1.0 |
| | Pump Pressure (bars): | 180 |
| | | |

REFERENCE SUBSTANCES:

| NO | Reference Substances/Related Substances | NQCL Code/Batch | Purity (%) |
|----|---|-----------------|------------|
| 1. | | | |
| 2. | | | |
| 3. | | | |
| 4. | | | |
| 5. | | | |