

CHROMATOGRAPHIC CONDITIONS:**ASSAY**

Column No: 100 Type of Column: waters Xterra C18
Column Temp (°C): 30
Detection λ (nm): 270 Injection Vol (μ L): 20

Mobile Phase: Composition (% v/v) & Ratios

Buffer : MeOH Flow Rate (mL/min): 1
50 : 50 Pump Pressure (bars): 244-250

DISSOLUTION

Column No: 100 Type of Column: waters Xterra C18
Column Temp (°C): 30
Detection λ (nm): 270 Injection Vol (μ L): 20

Mobile Phase: Composition (% v/v) & Ratios

Buffer : MeOH Flow Rate (mL/min): 1
50 : 50 Pump Pressure (bars): 244-250

REFERENCE SUBSTANCES:

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.	Lamivudine	22117702	100
2.	Zidovudine	WRS-Z01/06	99.4
3.	Nevirapine	C5028-15-195	98.8
4.			
5.			