CHROMATOGRAPHIC CONDITIONS:

		<u>ASSAY</u>		
Column No:	226	Type of Column:	Symmetry C18, 25cm	
Column Temp (°C):	40			
Detection λ (nm):	255	Injection Vol (μL):	10	
Mobile Phase: Compos	ition (% v/v) & Ratios		
BUFFER:ACN GRADIENT			Flow Rate (mL/min):	1.5
			Pump Pressure (bars):	224-60
		DISSOLUTI	<u>ON</u>	
Column No:	234	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	40	<u></u>		
Detection λ (nm):	255	Injection Vol (μL):	10	
Mobile Phase: Compos	ition (% v/v) & Ratios		
BUFFER:ACN_GRADIENT			Flow Rate (mL/min):	1.5
			Pump Pressure (bars):	224-60
REFERENCE SUBSTA	ANCES:			
210	(5.1	. 10.1	T C 1 /D : 1	0()

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			