CHROMATOGRAPHIC CONDITIONS:

		<u>ASSAY</u>		
Column No:	96	Type of Column:		
Column Temp (°C):	30	<u></u>		
Detection λ (nm):	254	Injection Vol (μL):	20	
Mobile Phase: Compos	ition (% v/v) & Ratios		
ACN TEA WATER 200 : 1 :700			Flow Rate (mL/min):	
200 . 1 .700			Pump Pressure (bars):	100-200
			<u></u>	
		DISSOLUTIO	<u>ON</u>	
Column No:	0	Type of Column:	N/A	
Column Temp (°C):	0	<u></u>		
Detection λ (nm):	N/A	Injection Vol (μL):	N/A	
Mobile Phase: Compos	ition (% v/v)) & Ratios	Flow Rate (mL/min):	N/A
n/a			Pump Pressure (bars):	N/A
				14//

REFERENCE SUBSTANCES:

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.	sulfamethoxazole	s12-6	99.2
2.	trimethoprim	t7-4	99.3
3.			
4.			
5.			