CHROMATOGRAPHIC CONDITIONS:

		<u>ASSAY</u>				
Column No:	226	Type of Column:	Symmetry C18, 25cm			
Column Temp (°C):	40					
Detection λ (nm):	255	Injection Vol (μL):	10			
Mobile Phase: Compos	ition (% v/v	7) & Ratios				
buffer : ACN			Flow Rate (mL/min):	1.5		
gradient			Pump Pressure (bars):	224-260		
<u>DISSOLUTION</u>						
Column No:	234	Type of Column:	Xterra RP 18, 25cm			
Column Temp (°C):	40					
Detection λ (nm):	255	Injection Vol (μL):	10			
Mobile Phase: Compos	ition (% v/v	y) & Ratios				
•	\ /	,	Flow Rate (mL/min):	1.5		
buffer : ACN gradient			Pump Pressure (bars):	224-260		
			<u> </u>			
REFERENCE SUBSTA	ANCES:					
		T				

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			