

DISSOLUTION

Standard Preparation for Dissolution:

Dissolution Conditions

	1 st Run	2 nd Run	3 rd Run
Dissolution Medium:	<u>BUFER pH 6.8</u>	<u>BUFER pH 6.8</u>	<u> </u>
Volume used:	<u>900</u>	<u>900</u>	<u> </u>
Apparatus:	<u>2</u>	<u>2</u>	<u> </u>
Rotations per minute:	<u>100</u>	<u>100</u>	<u> </u>
Time (min)	<u>45, 0, 0, 0, 0</u>	<u>45, 0, 0, 0, 0</u>	<u> </u>

Describe below any subsequent dilutions after the dissolution: