CHROMATOGRAPHIC CONDITIONS:

| | | ASSAY | | | |
|----------------------------|----------------|---------------------|-----------------------|---------|--|
| Column No: | 172 | Type of Column: | Symmetry C18, 25cm | | |
| Column Temp (°C): | 30 | | | | |
| Detection λ (nm): | 270 | Injection Vol (μL): | 20 | | |
| Mobile Phase: Compos | sition (% v/v) |) & Ratios | | | |
| Metnanol: Buffer ; 500:500 | | | Flow Rate (mL/min): | 1 | |
| | | | Pump Pressure (bars): | 285-290 | |
| | | | | | |
| | | | | | |
| | | DISSOLUTION | <u>ON</u> | | |
| Column No: | 172 | Type of Column: | Symmetry C18, 25cm | | |
| Column Temp (°C): | 30 | <u> </u> | | | |
| Detection λ (nm): | 270 | Injection Vol (μL): | 20 | | |
| Mobile Phase: Compos | sition (% v/v) | Nr Ratios | | | |
| • | , , , |) & Ratios | Flow Rate (mL/min): | | |
| Metnanol: Buffer ; 500:500 | | | | 1 | |
| | | | Pump Pressure (bars): | 285-290 | |
| | | | <u></u> | | |
| | | | | | |
| REFERENCE SUBST | ANCES: | | | | |
| | | | | | |

| NO | Reference Substances/Related Substances | NQCL Code/Batch | Purity (%) |
|----|---|-----------------|------------|
| 1. | | | |
| 2. | | | |
| 3. | | | |
| 4. | | | |
| 5. | | | |