CHROMATOGRAPHIC CONDITIONS:

		ASSA	<u>Y</u>	
Column No:	157	Type of Column:	Hyperclone BDS C18, 15cm	
Column Temp (°C):	30			
Detection λ (nm):	260	Injection Vol (μL)):20	
Mobile Phase: Composition (% v/v) & Ratios SODIUM ACETATE BUFFER pH 4.5			Flow Rate (mL/min): Pump Pressure (bars):	1 122-125
Column No.	457	DISSOLU'		
Column No:	157	_ Type of Column:	Hyperclone BDS C18, 15cm	
Column Temp (°C):	30	<u> </u>		
Detection λ (nm):	260	_ Injection Vol (μL)):20	
Mobile Phase: Compos	, ,	& Ratios	Flow Rate (mL/min):	1
		Pump Pressure (bars):	122-125	
REFERENCE SUBST	'ANCES:			
NO Defense Cul	-t/D-1-t	- 1 C. 11 NT	OCI Cada/Batala Dunitar	(0/)

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			