

**CHROMATOGRAPHIC CONDITIONS:****ASSAY**

Column No:	<u>169</u>	Type of Column:	<u>169</u>
Column Temp (°C):	<u>30</u>		
Detection λ (nm):	<u>254</u>	Injection Vol (μL):	<u>10</u>

Mobile Phase: Composition (% v/v) &amp; Ratios

<u>Buffer : Methanol (95:5)</u>	Flow Rate (mL/min):	<u>1.5</u>
	Pump Pressure (bars):	<u>150-170</u>

**DISSOLUTION**

Column No:	<u>169</u>	Type of Column:	<u>169</u>
Column Temp (°C):	<u>30</u>		
Detection λ (nm):	<u>254</u>	Injection Vol (μL):	<u>10</u>

Mobile Phase: Composition (% v/v) &amp; Ratios

<u>Buffer : Methanol (95:5)</u>	Flow Rate (mL/min):	<u>1.5</u>
	Pump Pressure (bars):	<u>150-170</u>

**REFERENCE SUBSTANCES:**

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.	Isoniazid	R013N0	1
2.			
3.			
4.			
5.			