CHROMATOGRAPHIC CONDITIONS:

		<u>ASSAY</u>		
Column No:	122	Type of Column:	Hyperclone BDS C18, 15cm	
Column Temp (°C):	30			
Detection λ (nm):	280	Injection Vol (μL):	10	
Mobile Phase: Compos	sition (% v/v	y) & Ratios		
Buffer:ACN - 30:70 pH 3.5			Flow Rate (mL/min):	1.0
pH 3.5			Pump Pressure (bars):	75-76
		DISSOLUTION	<u>ON</u>	
Column No:	122	Type of Column:	Hyperclone BDS C18, 15cm	
Column Temp (°C):	30			
Detection λ (nm):	280	Injection Vol (μL):	10	
Mobile Phase: Compos	sition (% v/v	y) & Ratios		
Buffer: ACN - 30:70 pH 3.5			Flow Rate (mL/min):	1.0
pH 3.5			Pump Pressure (bars):	75-76
REFERENCE SUBST	ANCES:			

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			