

**CHROMATOGRAPHIC CONDITIONS:****ASSAY**

Column No:	122	Type of Column:	Hyperclone BDS C18, 15cm
Column Temp (°C):	30		
Detection λ (nm):	214	Injection Vol (μL):	10

Mobile Phase: Composition (% v/v) &amp; Ratios

ACN:WATER	Flow Rate (mL/min):	1.0
23:77	Pump Pressure (bars):	65-75

**DISSOLUTION**

Column No:	122	Type of Column:	Hyperclone BDS C18, 15cm
Column Temp (°C):	30		
Detection λ (nm):	214	Injection Vol (μL):	10

Mobile Phase: Composition (% v/v) &amp; Ratios

ACN:WATER	Flow Rate (mL/min):	1.0
23:77	Pump Pressure (bars):	65-75

**REFERENCE SUBSTANCES:**

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			