

CHROMATOGRAPHIC CONDITIONS:

ASSAY

Column No: 179 Type of Column: Hyperclone BDS C18, 15cm
Column Temp (°C): 30
Detection λ (nm): 242 Injection Vol (μ L): 50

Mobile Phase: Composition (% v/v) & Ratios

acn-phosphate buffer gradient Flow Rate (mL/min): 1.2
Pump Pressure (bars): 185-186

DISSOLUTION

Column No: 179 Type of Column: Hyperclone BDS C18, 15cm
Column Temp (°C): 30
Detection λ (nm): 242 Injection Vol (μ L): 50

Mobile Phase: Composition (% v/v) & Ratios

acn-phosphate buffer gradient Flow Rate (mL/min): 1.2
Pump Pressure (bars): 185-186

REFERENCE SUBSTANCES:

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			