## **CHROMATOGRAPHIC CONDITIONS:**

		ASSAY		
Column No:	100	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	40			
Detection $\lambda$ (nm):	254	Injection Vol (μL):	20	_
Mobile Phase: Compos Water:ACN:TEA - 79.9:20.0	·	) & Ratios	Flow Rate (mL/min): Pump Pressure (bars):	1 120-121
		DISSOLUTIO	<u>ON</u>	
Column No:	100	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	40	<u></u>		
Detection $\lambda$ (nm):	254	Injection Vol (μL):	20	
Mobile Phase: Composition (% v/v) & Ratios  Water: ACN: TFA - 79.9:20.0:0.1			Flow Rate (mL/min): Pump Pressure (bars):	1 120-121
REFERENCE SUBST	ANCES:			

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			