## **CHROMATOGRAPHIC CONDITIONS:**

		<u>ASSAY</u>		
Column No:	168	Type of Column:	Symmetry C18, 25cm	
Column Temp (°C):	30	<u></u>		
Detection $\lambda$ (nm):	240	Injection Vol (μL):	20	
Mobile Phase: Compos	sition (% v/v)	& Ratios		
Methanol: 105g/L H3PO4 (80:20 V/V)			Flow Rate (mL/min):	1.0
	,		Pump Pressure (bars):	99-101
		DISSOLUTIO	<u>ON</u>	
Column No:	168	Type of Column:	Symmetry C18, 25cm	
Column Temp (°C):	30			
Detection $\lambda$ (nm):	240	Injection Vol (μL):	20	
14 1 D	(0/ / )	4 D		
Mobile Phase: Compos	sition (% v/v)	& Katios		
Methanol: 105g/l H3PO4 (80:20 V/V)			Flow Rate (mL/min):	1.0
			Pump Pressure (bars):	99-101
DEEEDENICE CHDCT	ANCEC.			
REFERENCE SUBST	ANCES:			

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			