

MICROBIAL COUNT

REFERENCE DOCUMENT: SOP NO. BIOL 007

MICROBIOLOGY LAB NO.	DATE RECEIVED	DATE TEST SET	DATE OF RESULTS
BIOL/001/2015	2015-12-17 06:23:55	23-Dec-2015	28-Dec-2015
SAMPLE PREPARATION			
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> 10g _____ 100ml 90 ml </div> <div style="text-align: center;">X</div> <div style="text-align: center;"> 10ml _____ 100ml 90 ml </div> <div style="text-align: center;">X</div> <div style="text-align: center;"> 1ml _____ 1 mlml Plating </div> <div style="margin-left: 20px;"> Replicates: 2 </div> </div>			
<div style="display: flex; justify-content: space-between;"> 0 0 </div>			
RESULTS			
		10 ¹ CFU 0	10 ² CFU 0
Nutrient Agar	Plate 1		
	Plate 2	0	0
	Average (A): CFU (Total Aerobic Microbial Count)	0	0
		0	0
Sabourauds Dextrose Agar	Plate 1	<100	0
	Plate 2		
	Average (B): CFU (Total Yeast Microbial Count)		
NB: Acceptance Criteria is interpreted as follows depending on route of administration – 10 ¹ cfu: maximum acceptable count = 20; 10 ² cfu: maximum acceptable count = 200; 10 ³ cfu: maximum acceptable count = 2000; and so forth.			
CONCLUSION: The Product		Complies	
		Does Not Comply	
Analyst:		Head, Biological Analysis Unit:	
Date:		Date:	
Analyst:		Signature:	
Date:			

TEST FOR SPECIFIED MICROORGANISMS

REFERENCE DOCUMENT: SOP NO. BIOL 007

BIOL/001/2015	2015-12-17 06:23:55	23-Dec-2015	28-Dec-2015
MICROBIOLOGY LAB NO.	DATE RECEIVED	DATE TEST SET	DATE OF RESULTS
SAMPLE PREPARATION			
<div>10g10ml1mlReplicates: 2</div> <div>100ml 90 ml100ml 90 ml1 mlml Plating</div>			
00			
RESULTS			
Microorganism	Test Media	Observation	Negative Control
		00	
		00	
		00	
		00	
		<1000	
Observation - Indicate whether there is growth/turbidity/colour change in the test media or Not.			
CONCLUSION: The Product	Complies	With the requirements of the Test for Specified Microorganisms.	
	Does Not Comply		
Analyst:		Head, Biological Analysis Unit:	
Date:		Date:	
Analyst:		Signature:	
Date:			