

MICROBIAL COUNT

REFERENCE DOCUMENT: SOP NO. BIOL 007

MICROBIOLOGY LAB NO.	DATE RECEIVED	DATE TEST SET	DATE OF RESULTS
SAMPLE PREPARATION			
RESULTS			
		10 ¹ CFU	10 ² CFU
Nutrient Agar	Plate 1		
	Plate 2		
	Average (A): CFU (Total Aerobic Microbial Count)		
Sabourauds Dextrose Agar	Plate 1		
	Plate 2		
	Average (B): CFU (Total Yeast Microbial Count)		
Negative Control			
Negative Control			
NB: Acceptance Criteria is interpreted as follows depending on route of administration – 10 ¹ cfu: maximum acceptable count = 20; 10 ² cfu: maximum acceptable count = 200; 10 ³ cfu: maximum acceptable count = 2000; and so forth.			
CONCLUSION: The Product		Complies	
		Does Not Comply	
Analyst:		Head, Biological Analysis Unit:	
Date:		Date:	
Analyst:		Signature:	
Date:			

TEST FOR SPECIFIED MICROORGANISMS
REFERENCE DOCUMENT: SOP NO. BIOL 007

MICROBIOLOGY LAB NO.	DATE RECEIVED	DATE TEST SET	DATE OF RESULTS
SAMPLE PREPARATION			
RESULTS			
Microorganism	Test Media	Observation	Negative Control
Observation – Indicate whether there is growth/turbidity/colour change in the test media or Not.			
CONCLUSION: The Product		Complies	With the requirements of the Test for Specified Microorganisms.
		Does Not Comply	
Analyst:		Head, Biological Analysis Unit:	
Date:		Date:	
Analyst:		Signature:	
Date:			