

CHROMATOGRAPHIC CONDITIONS:**ASSAY**

Column No:	<u>163</u>	Type of Column:	<u>Hypersil BDS Cyano, 25cm</u>
Column Temp (°C):	<u>30</u>		
Detection λ (nm):	<u>265</u>	Injection Vol (μL):	<u>20</u>

Mobile Phase: Composition (% v/v) & Ratios

<u>BUFFER OF 1ML H3PO4 IN 1L H2O:ACN</u>	Flow Rate (mL/min):	<u>1.8</u>
<u>68:32</u>	Pump Pressure (bars):	<u>122-125</u>

DISSOLUTION

Column No:	<u>163</u>	Type of Column:	<u>Hypersil BDS Cyano, 25cm</u>
Column Temp (°C):	<u>30</u>		
Detection λ (nm):	<u>265</u>	Injection Vol (μL):	<u>20</u>

Mobile Phase: Composition (% v/v) & Ratios

<u>BUFFER OF 1ML H3PO4 IN 1L H2O:ACN</u>	Flow Rate (mL/min):	<u>1.8</u>
<u>68:32</u>	Pump Pressure (bars):	<u>122-125</u>

REFERENCE SUBSTANCES:

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			