

CHROMATOGRAPHIC CONDITIONS:**ASSAY**

Column No:	155	Type of Column:	Eclipse plus C8, 4.6 * 250mm
Column Temp (°C):	25		
Detection λ (nm):	210	Injection Vol (μL):	10

Mobile Phase: Composition (% v/v) & Ratios

BUFFER: SOLVENT MIXTURE	Flow Rate (mL/min):	1.5
450:550	Pump Pressure (bars):	215-225

DISSOLUTION

Column No:	155	Type of Column:	Eclipse plus C8, 4.6 * 250mm
Column Temp (°C):	25		
Detection λ (nm):	210	Injection Vol (μL):	10

Mobile Phase: Composition (% v/v) & Ratios

BUFFER: SOLVENT MIXTURE	Flow Rate (mL/min):	1.5
450:550	Pump Pressure (bars):	215-225

REFERENCE SUBSTANCES:

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.	LOPINAVIR	L20-2	99.8
2.	RITONAVIR	R14-2	99.4
3.			
4.			
5.			