


M19 -7

Certificate of Analysis

ISO GUIDE 34
ACCLASS Cert# AR-1470

ISO/IEC 17025
ACCLASS Cert# AT-1467

METFORMIN HYDROCHLORIDE CERTIFIED REFERENCE MATERIAL



CERTIFIED PURITY: 99.8%, $U_{\text{CRM}} = \pm 0.4\%$ $k = 2.18$
(Mass Balance/as is basis)

NOMINAL PACKAGE SIZE: 500mg

CATALOG #: PHR1084

LOT #: LRAA1730

CERTIFICATE VERSION: LRAA1730.1

ISSUE DATE: 17 January 2014

*Note: Certificates may be updated due to Pharmacopeial Lot changes or the availability of new data.
Check our website at: www.sigma-aldrich.com for the most current version.*

CRM EXPIRATION: 12 Months from Receipt (Proper Storage and Handling Required).

RECEIPT DATE: _____

Note: this space is provided for convenience only and its use is not required.

STORAGE: Store in a Refrigerator, keep container tightly closed. Attachment of a 20 mm aluminum crimp seal recommended for unused portions.

CHEMICAL FORMULA: $\text{C}_4\text{H}_{12}\text{ClN}_5$

MW: 165.6


PHYSICAL DESCRIPTION: White powder in amber vial

CAS #: 1115-70-4

HAZARDS: Read Safety Data Sheet before using. All chemical reference materials should be considered potentially hazardous and should be used only by qualified laboratory personnel.



SIGMA-ALDRICH®



INSTRUCTIONS FOR USE: Do not dry, use on the as is basis. The internal pressure of the container may be slightly different from the atmospheric pressure at the user's location. Open slowly and carefully to avoid dispersion of the material. This material is intended for R&D use only. Not for drug, household or other uses.

TRACEABILITY ASSAY

Comparative assay demonstrates direct traceability to Pharmacopeial Standards
Specification: 98.5 to 101.0% (USP)

METHOD: HPLC (ref.: Glyburide and Metformin Hydrochloride Tablets; USP36)

Column: Ascentis C18, 4.6 x 250mm, 5µm

Mobile Phase: 0.5g/L Sodium-1-heptanesulfonate and 0.5g/L NaCl in Water, Acetonitrile (90:10)

Flow: 1mL/min

Column Temperature: 30°C

Injection: 3.5µl

Detector: 218nm

ASSAY vs. USP REFERENCE STANDARD (as is basis)

ASSAY VALUE

99.98%

vs. USP LOT

J0L465

Labeled Content = 0.998mg/mg

ASSAY vs. EP CRS (as is basis)

ASSAY VALUE

99.8%

vs. EP BATCH

3.1

Labeled Content = None

Assigned Content = 99.5%*

ASSAY vs. BP CRS (as is basis)

ASSAY VALUE

99.5%


vs. BP BATCH

3267

Labeled Content = 99.8%

*The assigned content of the EP CRS was determined by assay against the USP Reference Standard and the BP CRS

Representative Chromatogram from Lot: LRAA1730 Analysis



PURITY DETERMINATION BY MASS BALANCE

CHROMATOGRAPHIC IMPURITY ANALYSIS

METHOD: HPLC (ref.: Metformin Hydrochloride; USP34)

Column: Supelcosil LC- SCX, 4.6 x 250mm, 5 μ m

Mobile Phase: Monobasic ammonium phosphate in water, pH 3.0 with phosphoric acid

Flow: 1mL/min

Column Temperature: 30°C

Injection: 2 μ l


Detector: 218nm

Impurities Detected:

Impurity 1:	0.001%
1-Cyanoguanidine:	0.002%
Impurity 2:	0.004%

Total Impurities: 0.006%

Representative Chromatogram from Lot: LRAA1730 Impurity Analysis



RESIDUAL SOLVENTS

Method: GC-MS Headspace (ref.: Residual Solvents <467>, USP34)

Column: DB-1301

Carrier gas: He

Flow: 1.2mL/min

Split Ratio: 1:5

Injection/Temperature: 1mL/250°C

Temperature Program: 40°C for 20min, 10°C/min to 240°C, hold 20min

Solvents Detected: None

LOSS ON DRYING/VOLATILES

Method: 105°C

Mean of three samples, Loss = 0.2%

RESIDUE ANALYSIS

Method: Sulfated Ash

Sample Size: ~1g

Mean of three determinations, Residue = 0.02%


CERTIFIED PURITY BY MASS BALANCE [100% - Impurities (normalized)]

99.8% $U_{\text{crm}} = \pm 0.4\%$, $k = 2.18$
(as is basis)

IDENTIFICATION TESTS

INFRARED SPECTROPHOTOMETRY (Comparative identification analysis demonstrates direct traceability to Pharmacopeial standards)

FTIR Comparison of Metformin Hydrochloride




— Metformin Hydrochloride USP Lot JOL465
— Metformin Hydrochloride EP Batch 3.1
— Metformin Hydrochloride BP Batch 3267
— Metformin Hydrochloride RTC Lot LRAA1730

PERKIN ELMER SPECTRUM 100
SAMPLING: UNIVERSAL ATR

¹H NMR (Data provided by an external laboratory; not in scope of accreditation)

LRAA1730 Metformin Hydrochloride



Consistent with structure

ELEMENTAL ANALYSIS (Data provided by an external laboratory; not in scope of accreditation)

Exeter Analytical 440 Elemental Analyzer

Combustion method

%	Theoretical	Result 1	Result 2	Mean
C	29.00	29.32	29.29	29.31
H	7.30	6.97	7.00	6.99
N	42.28	42.28	42.36	42.32

MELTING RANGE

Specification: 222°C to 226°C (EP)

Mettler Toledo FP900 Thermosystem with FP81 Measuring Cell

Mean of three measurements = 223.1° – 225.5°C

HOMOGENEITY ASSESSMENT

Homogeneity was assessed in accordance with ISO Guide 35. Completed units were sampled using a random stratified sampling protocol. The results of chemical analysis were then compared by Single Factor Analysis of Variance (ANOVA). The uncertainty due to homogeneity was derived from the ANOVA. Heterogeneity was not detected under the conditions of the ANOVA.

Analytical Method: HPLC

Sample size: ~25mg

UNCERTAINTY STATEMENT


Uncertainty values in this document are expressed as Expanded Uncertainty (U_{crm}) corresponding to the 95% confidence interval. U_{crm} is derived from the combined standard uncertainty multiplied by the coverage factor k , which is obtained from a t -distribution and degrees of freedom. The components of combined standard uncertainty include the uncertainties due to characterization, homogeneity, long term stability, and short term stability (transport). The components due to stability are generally considered to be negligible unless otherwise indicated by stability studies.

STABILITY ASSESSMENT

Significance of the stability assessment will be demonstrated if the analytical result of the study and the range of values represented by the Expanded Uncertainty do not overlap the result of the original assay and the range of its values represented by the Expanded Uncertainty. The method employed will usually be the same method used to characterize the assay value in the initial evaluation.

Long Term Stability Evaluation - An assessment, or re-test, versus a Compendial Reference Standard may be scheduled, within the 3 year anniversary date of a release of a Secondary Standard. The re-test interval will be determined on a case-by-case basis.

Short Term Stability Study - It is useful to assess stability under reasonably anticipated, short term transport conditions by simulating exposure of the product to humidity and temperature stress. This type of study is conducted under controlled conditions of elevated temperature and humidity.



Operations Manager



QA Supervisor

APPENDIX

Original Release Date: 17 January 2014

Manufactured and certified by Sigma-Aldrich RTC, Inc.
2931 Soldier Springs Rd, Laramie WY, USA 82070
(Phone): 1-307-742-5452 (Fax): 1-855-831-9212
email: RTCTechGroup@sial.com

