

DISSOLUTION

Standard Preparation for Dissolution:

Dissolution Conditions

	1 st Run	2 nd Run	3 rd Run
Dissolution Medium:	SODIUM PHOSPHATE BUFFER pH 6.8	SODIUM PHOSPHATE BUFFER pH 6.8	SODIUM PHOSPHATE BUFFER pH 6.8
Volume used:	900	900	
Apparatus:	2	2	
Rotations per minute:	100	100	
Time (min)	45, 0, 0, 0, 0	45, 0, 0, 0, 0	

Describe below any subsequent dilutions after the dissolution: