## **CHROMATOGRAPHIC CONDITIONS:**

		<u>ASSAY</u>			
Column No:	103	Type of Column:	Hyperclone BDS C18, 15cm		
Column Temp (°C):	30				
Detection $\lambda$ (nm):	275	Injection Vol (μL):	20		
Mobile Phase: Compos <u>SOL A: SOL B</u> GRADIENT	ition (% v/v)	& Ratios	Flow Rate (mL/min): Pump Pressure (bars):	2 162-164	
DISSOLUTION  Column No: 103 Type of Column: Hyperclone BDS C18, 15cm					
Column Temp (°C):	30	<u> </u>			
Detection $\lambda$ (nm):	237	Injection Vol (μL):	50		
Mobile Phase: Compos H20:ACN:TFA B5:35.0 PH 3.0	ition (% v/v)	& Ratios	Flow Rate (mL/min): Pump Pressure (bars):	1.5 156-157	
REFERENCE SUBSTA	ANCES:				

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			