## **CHROMATOGRAPHIC CONDITIONS:**

		<u>ASSAY</u>		
Column No:	234	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	35			
Detection λ (nm):	254	Injection Vol (μL):	20	
Mobile Phase: Composi	•		Flow Rate (mL/min):	1.8
H2O:ACN: TEA, 1400:400:2 adjusted to pH 5.9 with GAA and made upto 2000ml with H2O			Pump Pressure (bars):	148-159
		DISSOLUTIO		
Column No:	234	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	35		,	
Detection λ (nm):	254	Injection Vol (μL):	20	
Mobile Phase: Composi	tion (% v/v	7) & Ratios		
H2O:ACN: TFA 1400:400:2 adjusted to pt. 5.9 with GAA and made upto 2000ml			Flow Rate (mL/min):	1.8
With H20			Pump Pressure (bars):	148-159
			<u></u>	

## **REFERENCE SUBSTANCES:**

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			