

CHROMATOGRAPHIC CONDITIONS:

ASSAY

Column No: 250 Type of Column: Thermo Betasil C1, 25cm
Column Temp (°C): 30
Detection λ (nm): 270 Injection Vol (μ L): 10

Mobile Phase: Composition (% v/v) & Ratios

BUFFER:MEOH Flow Rate (mL/min): 1
500:500 Pump Pressure (bars): 195-205

DISSOLUTION

Column No: 250 Type of Column: Thermo Betasil C1, 25cm
Column Temp (°C): 30
Detection λ (nm): 270 Injection Vol (μ L): 10

Mobile Phase: Composition (% v/v) & Ratios

BUFFER:MEOH Flow Rate (mL/min): 1
500:500 Pump Pressure (bars): 195-205

REFERENCE SUBSTANCES:

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			