

CHROMATOGRAPHIC CONDITIONS:**ASSAY**

Column No: 132 Type of Column: Hypersil BDS Cyano, 15 cm
Column Temp (°C): 35
Detection λ (nm): 254 Injection Vol (μ L): 20

Mobile Phase: Composition (% v/v) & Ratios

soln B: soln C Flow Rate (mL/min): 1.5
Pump Pressure (bars): 70-75

DISSOLUTION

Column No: _____ Type of Column: _____,
Column Temp (°C): _____
Detection λ (nm): _____ Injection Vol (μ L): _____

Mobile Phase: Composition (% v/v) & Ratios

Flow Rate (mL/min): _____
Pump Pressure (bars): _____

REFERENCE SUBSTANCES:

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			