

**CHROMATOGRAPHIC CONDITIONS:****ASSAY**

Column No:	<u>97</u>	Type of Column:	<u>c18,25*4.6mm,5um</u>
Column Temp (°C):	<u>30</u>		
Detection λ (nm):	<u>254</u>	Injection Vol (μL):	<u>10</u>

Mobile Phase: Composition (% v/v) &amp; Ratios

<u>ACN:GAA:WATER(12:1:87)</u>	Flow Rate (mL/min):	<u>2</u>
	Pump Pressure (bars):	<u>150-180</u>

**DISSOLUTION**

Column No:	<u>0</u>	Type of Column:	<u>N/A</u>
Column Temp (°C):	<u>0</u>		
Detection λ (nm):	<u>N/A</u>	Injection Vol (μL):	<u>N/A</u>

Mobile Phase: Composition (% v/v) &amp; Ratios

<u>n/a</u>	Flow Rate (mL/min):	<u>N/A</u>
	Pump Pressure (bars):	<u>N/A</u>

**REFERENCE SUBSTANCES:**

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.	Trimethoprim	VWS/135/15	99.3
2.	Sulfadiazine	VWS/120/15	99.5
3.			
4.			
5.			