

**CHROMATOGRAPHIC CONDITIONS:****ASSAY**

Column No:	<u>125</u>	Type of Column:	<u>Hyperclone BDS C18, 15cm</u>
Column Temp (°C):	<u>35</u>		
Detection λ (nm):	<u>261</u>	Injection Vol (μL):	<u>20</u>

Mobile Phase: Composition (% v/v) &amp; Ratios

<u>MeOH:ACN: bUFFER, 20:10:70 PH 5.0</u>	Flow Rate (mL/min):	<u>1</u>
	Pump Pressure (bars):	<u>74-80</u>

**DISSOLUTION**

Column No:	<u>125</u>	Type of Column:	<u>Hyperclone BDS C18, 15cm</u>
Column Temp (°C):	<u>35</u>		
Detection λ (nm):	<u>261</u>	Injection Vol (μL):	<u>20</u>

Mobile Phase: Composition (% v/v) &amp; Ratios

<u>MeOH:ACN: bUFFER, 20:10:70 PH 5.0</u>	Flow Rate (mL/min):	<u>1</u>
	Pump Pressure (bars):	<u>74-80</u>

**REFERENCE SUBSTANCES:**

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			