CHROMATOGRAPHIC CONDITIONS:

		ASSAY				
Column No:	256	Type of Column:	C18			
Column Temp (°C):	30	<u></u>				
Detection λ (nm):	254	Injection Vol (μL):	10			
Mobile Phase: Compos	sition (% v/v	y) & Ratios				
WATER:ACN:TEA =			Flow Rate (mL/min):	1.0		
1400:400:2 in 2000 ML, pH 5.9			Pump Pressure (bars):	145-150		
			<u></u>			
		DISCOLUTIV	av.			
<u>DISSOLUTION</u>						
Column No:	256	Type of Column:	C18			
Column Temp (°C):	30					
Detection λ (nm):	254	Injection Vol (μL):	10			
Mobile Phase: Compos	sition (% v/v	y) & Ratios				
WATER·ACN·TEA =	·	Flow Rate (mL/min):	1.0			
1400:400:2 in 2000 ML, pF	1 5.9	Pump Pressure (bars):	145-150			

REFERENCE SUBSTANCES:

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.	SULFAMETHOXAZOLE	SMX15082073	99.02
2.	TRIMETHOPRIM	VWS/135/15	99
3.			
4.			
5.			