## **CHROMATOGRAPHIC CONDITIONS:**

		<u>ASSAY</u>		
Column No:	100	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	30			
Detection $\lambda$ (nm):	260	Injection Vol (μL):	20	
Mobile Phase: Compos	sition (% v/v)	& Ratios	Elavy Data (m. I. / main)	1
acn:buffer(40:60)			Flow Rate (mL/min):	1
			Pump Pressure (bars):	200-220
			<u></u>	
		DISSOLUTI	<u>ON</u>	
Column No:	100	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	30			
Detection $\lambda$ (nm):	260	Injection Vol (μL):	20	
M.1.1. Dl C		A Dada		
Mobile Phase: Compos	31tion (% V/V)	) & Kaπos	Flow Rate (mL/min):	
acn:buffer(40:60)				
			Pump Pressure (bars):	200-220
REFERENCE SUBST	ANCES:			
210 7 6 6 1	. /5 :	. 10.1 .		0/)

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			