

**DISSOLUTION**

Standard Preparation for Dissolution:

**Dissolution Conditions**

|                       | 1 <sup>st</sup> Run | 2 <sup>nd</sup> Run | 3 <sup>rd</sup> Run |
|-----------------------|---------------------|---------------------|---------------------|
| Dissolution Medium:   | 0.1MHCl             | 0.1MHCl             |                     |
| Volume used:          | 900                 | 900                 |                     |
| Apparatus:            | 2                   | 2                   |                     |
| Rotations per minute: | 50                  | 50                  |                     |
| Time (min)            | 45, 0, 0, 0, 0      | 45, 0, 0, 0, 0      |                     |

Describe below any subsequent dilutions after the dissolution: