CHROMATOGRAPHIC CONDITIONS:

		<u>ASSAY</u>		
Column No:	171	Type of Column:	Symmetry C18, 25cm	
Column Temp (°C):	30			
Detection λ (nm):	270	Injection Vol (μL):	20	
Mobile Phase: Compos	, -	r) & Ratios	Flow Rate (mL/min):	1
INICOLLAGE LATE BOLLER 30:30			Pump Pressure (bars): 210-230	
			<u> </u>	
		DISSOLUTION	<u>ON</u>	
Column No:	171	Type of Column:	Symmetry C18, 25cm	
Column Temp (°C):	30	<u></u>		
Detection λ (nm):	270	Injection Vol (μL):	20	
Mobile Phase: Compos	sition (% v/v	r) & Ratios		
MEOH:ACETATE BUFFER 50:50			Flow Rate (mL/min):	1
			Pump Pressure (bars):	210-230
REFERENCE SUBST	ANCES:			
NO Reference Sub	stances/Rela	ated Substances NO	CL Code/Batch Purity (%)

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			