CHROMATOGRAPHIC CONDITIONS:

		<u>ASSAY</u>		
Column No:	233	Type of Column:	3.9MM*30CM	
Column Temp (°C):	30			
Detection λ (nm):	254	Injection Vol (μL):	20	
Mobile Phase: Compos	sition (% v/v)	& Ratios		
ACN TEA WATER			Flow Rate (mL/min):	2
200 : 1 :700			Pump Pressure (bars):	100-200
			<u> </u>	
		DISSOLUTIO	ON	
Column No:	258	Type of Column:	4.6*250,C18	
Column Temp (°C):	30			
Detection λ (nm):	254	Injection Vol (μL):	20	
Mobile Phase: Compos	sition (% v/v)	& Ratios		
ACN WATER TEA			Flow Rate (mL/min):	1.0
200 :700 :1			Pump Pressure (bars):	100-200

REFERENCE SUBSTANCES:

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.	sulfamethoxazole	s12-6	99.2
2.	trimethoprim	t7-4	99.3
3.			
4.			
5.			