## **CHROMATOGRAPHIC CONDITIONS:**

5.

		<u>ASSA</u>	<u>\Y</u>	
Column No:	196	Type of Column	: Xterra RP 18, 25cm	
Column Temp (°C):	30			
Detection $\lambda$ (nm):	214	Injection Vol (μΙ	L):20	
Mobile Phase: Compos	ition (% v/	v) & Ratios		
ACN:WATER(23:77)			Flow Rate (mL/min):	1
			Pump Pressure (bars):	200-230
		DISSOLU		
Column No:	193	Type of Column	: Xterra RP 18, 25cm	
Column Temp (°C):	30			
Detection $\lambda$ (nm):	214	Injection Vol (μΙ	L): <u>20</u>	
Mobile Phase: Compos	ition (% v/v	v) & Ratios		
ACN:WATER(23:77)	•	Flow Rate (mL/min):	1	
- ( - )			Pump Pressure (bars):	163-200
REFERENCE SUBSTA	ANCES:			
NO Reference Sub	stances / Rel	ated Substances N	JOCL Code/Batch Purity (	(%)

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			