

CHROMATOGRAPHIC CONDITIONS:

ASSAY

| | | | |
|---------------------------|-----|---------------------------|--------------------|
| Column No: | 226 | Type of Column: | Symmetry C18, 25cm |
| Column Temp (°C): | 40 | | |
| Detection λ (nm): | 255 | Injection Vol (μ L): | 10 |

Mobile Phase: Composition (% v/v) & Ratios

| | | |
|--------------------------|-----------------------|--------|
| buffer : ACN gradient | Flow Rate (mL/min): | 1.5 |
| | Pump Pressure (bars): | 224-60 |

DISSOLUTION

| | | | |
|---------------------------|-----|---------------------------|--------------------|
| Column No: | 234 | Type of Column: | Xterra RP 18, 25cm |
| Column Temp (°C): | 40 | | |
| Detection λ (nm): | 255 | Injection Vol (μ L): | 10 |

Mobile Phase: Composition (% v/v) & Ratios

| | | |
|--------------------------|-----------------------|--------|
| buffer : ACN gradient | Flow Rate (mL/min): | 1.5 |
| | Pump Pressure (bars): | 224-60 |

REFERENCE SUBSTANCES:

| NO | Reference Substances/Related Substances | NQCL Code/Batch | Purity (%) |
|----|---|-----------------|------------|
| 1. | | | |
| 2. | | | |
| 3. | | | |
| 4. | | | |
| 5. | | | |