

CHROMATOGRAPHIC CONDITIONS:

ASSAY

Column No:	165	Type of Column:	Hypersil BDS Cyano, 25cm
Column Temp (°C):	40		
Detection λ (nm):	200	Injection Vol (μL):	25

Mobile Phase: Composition (% v/v) & Ratios

	Flow Rate (mL/min):	1
	Pump Pressure (bars):	90-100

DISSOLUTION

Column No:	165	Type of Column:	Hypersil BDS Cyano, 25cm
Column Temp (°C):	40		
Detection λ (nm):	200	Injection Vol (μL):	25

Mobile Phase: Composition (% v/v) & Ratios

Buffer: ACN	Flow Rate (mL/min):	1
50 :50	Pump Pressure (bars):	90-100

REFERENCE SUBSTANCES:

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			