CHROMATOGRAPHIC CONDITIONS:

		<u>ASSAY</u>		
Column No:	191	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	30			
Detection λ (nm):	260	Injection Vol (μL):	20	
Mobile Phase: Compos	sition (% v/v	y) & Ratios		
Buffer (0.05M NaH2PO4):ACN, 60:40			Flow Rate (mL/min):	1
			Pump Pressure (bars):	157-162
		DISSOLUT	<u>ION</u>	
Column No:	191	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	30			
Detection λ (nm):	260	Injection Vol (μL):	20	
Mohile Phase: Compos	sition (% y/x	v) & Ratios		
Mobile Phase: Composition (% v/v) & Ratios			Flow Rate (mL/min):	1
Buffer (0.05M NaH2PO4):ACN, 60:40			Pump Pressure (bars):	157-162
REFERENCE SUBST.	ANCES:			
NO Reference Sub	ctancos /Rol	ated Substances NO	CI Codo/Batch Purity (%)

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			