

CHROMATOGRAPHIC CONDITIONS:

ASSAY

Column No: 100 Type of Column: Xterra RP 18, 25cm
Column Temp (°C): 30
Detection λ (nm): 220 Injection Vol (μL): 20

Mobile Phase: Composition (% v/v) & Ratios

MEOH:BUFFER Flow Rate (mL/min): 1.5
1: 19 Pump Pressure (bars): 200-220

DISSOLUTION

Column No: 100 Type of Column: Xterra RP 18, 25cm
Column Temp (°C): 30
Detection λ (nm): 220 Injection Vol (μL): 20

Mobile Phase: Composition (% v/v) & Ratios

meoh:buffer Flow Rate (mL/min): 1.5
1: 19 Pump Pressure (bars): 200-220

REFERENCE SUBSTANCES:

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			