## **CHROMATOGRAPHIC CONDITIONS:**

		<u>ASSAY</u>		
Column No:	100	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	40	<u>_</u>		
Detection $\lambda$ (nm):	294,310	_ Injection Vol (μL):	20	
Mobile Phase: Composi KH2PO4 BUFFER:ACN 70:30	ition (% v/v)	& Ratios	Flow Rate (mL/min): Pump Pressure (bars):	1.2 145-150
		DISSOLUTIO	<u> </u>	
Column No:	100	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	40	<u>_</u>		
Detection $\lambda$ (nm):	294,310	_ Injection Vol (μL):	20	
Mobile Phase: Composition (% v/v) & Ratios KH2PO4 BUFFER:ACN 70:30			Flow Rate (mL/min): Pump Pressure (bars):	1.2 145-150
REFERENCE SUBSTA	ANCES:			

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			