CHROMATOGRAPHIC CONDITIONS:

		<u>ASSAY</u>		
Column No:	94	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	30			
Detection λ (nm):	261	Injection Vol (μL):	20	
Mobile Phase: Compos buffer:acn(60:40)	sition (% v/v	v) & Ratios	Flow Rate (mL/min): Pump Pressure (bars):	1 150-160
		DISSOLUT	<u>ION</u>	
Column No:	94	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	30			
Detection λ (nm):	261	Injection Vol (μL):	20	
Mobile Phase: Compos	sition (% v/v	v) & Ratios		
buffer:methanol(60:40)			Flow Rate (mL/min):	1
			Pump Pressure (bars):	130-160
			<u></u>	
REFERENCE SUBST	ANCES:			
NO Reference Sub	stances/Rel	ated Substances NO	CL Code/Batch Purity (%)

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			