

**CHROMATOGRAPHIC CONDITIONS:**

**ASSAY**

Column No:	<u>123</u>	Type of Column:	<u>Hyperclone BDS C18, 15cm</u>
Column Temp (°C):	<u>30</u>		
Detection λ (nm):	<u>210,295</u>	Injection Vol (μL):	<u>10</u>

Mobile Phase: Composition (% v/v) & Ratios

<u>pH 3.5 buffer:ACN:MeOH:THF</u>	Flow Rate (mL/min):	<u>2.2</u>
<u>48:42:5:5</u>	Pump Pressure (bars):	<u>106-10</u>

**DISSOLUTION**

Column No:	<u>123</u>	Type of Column:	<u>Hyperclone BDS C18, 15cm</u>
Column Temp (°C):	<u>30</u>		
Detection λ (nm):	<u>210,295</u>	Injection Vol (μL):	<u>10</u>

Mobile Phase: Composition (% v/v) & Ratios

<u>pH 3.5 buffer:ACN:MeOH:THF</u>	Flow Rate (mL/min):	<u>2.2</u>
<u>48:42:5:5</u>	Pump Pressure (bars):	<u>106-10</u>

**REFERENCE SUBSTANCES:**

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			