## **CHROMATOGRAPHIC CONDITIONS:**

		<u>ASSAY</u>		
Column No:	226	Type of Column:	Symmetry C18, 25cm	
Column Temp (°C):	40			
Detection $\lambda$ (nm):	1.5	Injection Vol (μL):	10	
Mobile Phase: Compos	sition (% v/v	v) & Ratios		
buffer : ACN			Flow Rate (mL/min):	1.5
gradient			Pump Pressure (bars):	224-260
		DISSOLUTIO	ON	
Column No:	234	Type of Column:	Xterra RP 18, 25cm	
Column Temp (°C):	40			
Detection $\lambda$ (nm):	1.5	Injection Vol (μL):	10	
Mobile Phase: Compos	sition (% v/v	7) & Ratios		
-	,	,	Flow Rate (mL/min):	1.5
buffer : ACN gradient			Pump Pressure (bars):	224-260
REFERENCE SUBST	ANCES:			

NO	Reference Substances/Related Substances	NQCL Code/Batch	Purity (%)
1.			
2.			
3.			
4.			
5.			